

Sizing Guide




This guide will assist in selecting the correct respirator size for your face. Begin with measuring the vertical line as shown to determine the dimension between the nose-bridge area to the chin.


Vertical Line Measurement



Once you determine the vertical dimension, measure the horizontal line. In order to obtain a proper fit and airtight seal, your measurements cannot be greater than the dimensions shown below. Refer to images below of interior of the Comfort-Air®NxMD respirators.



400NxMD Small / Medium
Vertical Line - 3.75"
Horizontal Line - 4.55"



400NxMD Medium / Large
Vertical Line - 4.5"
Horizontal Line - 3.75"

Scan to Buy

